
















Radius PPG™

Wireless SpO2 Sensors

Images	2-3
en English	4-7
fr French	8-11
de German	12-15
it Italian	16-19
es Spanish	20-23
sv Swedish	24-27
nl Dutch	28-31
da Danish	32-35
pt Portuguese	36-39
zh Chinese	40-42
ja Japanese	43-46
fi Finnish	47-50
no Norwegian	51-54
cs Czech	55-58
hu Hungarian	59-62
pl Polish	63-66
ro Romanian	67-70
sk Slovak	71-74
tr Turkish	75-78
el Greek	79-82
ru Russian	83-86
ko Korean	87-90
ar Arabic	94-91
Performance Specifications	95


Fig. 8a


Fig. 8b


Fig. 9


Fig. 10


Fig. 11


Fig. 12

Radius PPG™

Wireless SpO2 Sensors

en

DIRECTIONS FOR USE

-  Single patient use only

Not made with natural rubber latex

Non-sterile

Prior to using this product the user should read and understand the Operator's Manual for the device/monitor, the Directions for Use for the Radius PPG chip and receiver, and this Directions for Use.

INDICATIONS

Masimo Radius PPG™ is intended for the non-invasive continuous monitoring of functional oxygen saturation of arterial hemoglobin (SpO₂), pulse rate (PR).

Masimo Radius PPG is indicated for the continuous monitoring of functional arterial oxygen saturation of hemoglobin (SpO₂) and pulse rate (PR) for use with adult, pediatric and neonatal patients during both no motion and motion conditions and for patients who are well or poorly perfused in hospital, hospital-type facilities and home environments.

Devices with Masimo technology are only to be used with Masimo sensors and cables.

CONTRAINDICATIONS

The Radius PPG is contraindicated for patients who exhibit allergic reactions to foam rubber products and/or adhesive tape.

DESCRIPTION

Radius PPG consists of three parts:

- **Radius PPG wireless receiver**
(shipped separately, see *Directions for Use for the Radius PPG chip and receiver*)
- **Radius PPG reusable chip**
(shipped separately, see *Directions for Use for the Radius PPG chip and receiver*)
- **Radius PPG adhesive sensor**

Radius PPG is a wireless sensor for use with devices containing Masimo technology MX Version 7.14.8.x. or higher. Consult individual device manufacturers for compatibility of particular device and sensor models. Each device manufacturer is responsible for determining whether their device is compatible with each sensor model.

WARNINGS, CAUTIONS AND NOTES

- All sensors and cables are designed for use with specific monitors. Verify the compatibility of the monitor, cable and sensor before use, otherwise degraded performance and/or patient injury can result.
- The sensor should be free of visible defects, discoloration and damage. If the sensor is discolored or damaged, discontinue use. Never use a damaged sensor or one with exposed electrical circuitry.
- The site must be checked frequently or per clinical protocol to ensure adequate adhesion, circulation, skin integrity and correct optical alignment.
- Do not leave the sensor components unattended around children. Small items may become choking hazards.
- Exercise caution with poorly perfused patients; skin erosion and pressure necrosis can be caused when the sensor is not frequently moved. Assess site as frequently as every (1) hour with poorly perfused patients and move the sensor if there are signs of tissue ischemia.
- Circulation distal to the sensor site should be checked routinely.
- The sensor attachment strap site must be checked frequently or per clinical protocol to ensure adequate securement, circulation and skin integrity.
- During low perfusion, the sensor site needs to be assessed frequently for signs of tissue ischemia, which can lead to pressure necrosis.
- With very low perfusion at the monitored site, the reading may read lower than core arterial oxygen saturation.
- Do not use tape to secure the sensor to the site; this can restrict blood flow and cause inaccurate readings. Use of additional tape can cause skin damage, and/or pressure necrosis or damage the sensor.
- Sensors applied too tightly or that become tight due to edema will cause inaccurate readings and can cause pressure necrosis.
- Misapplied sensors or sensors that become partially dislodged may cause incorrect measurements.
- Misapplications due to wrong sensor types can cause inaccurate or no readings.
- Venous congestion may cause under reading of actual arterial oxygen saturation. Therefore, assure proper venous outflow from monitored site. Sensor should not be below heart level (e.g. sensor on hand of a patient in a bed with arm dangling to the floor).
- Venous pulsations may cause erroneous low SpO₂ readings (e.g. tricuspid valve regurgitation).
- The pulsations from intra-aortic balloon support can affect the pulse rate displayed on the oximeter. Verify patient's pulse rate against the ECG heart rate.
- Carefully route the sensor cable reduce the possibility of patient entanglement or strangulation.
- Avoid placing the sensor on any extremity with an arterial catheter or blood pressure cuff.
- If using pulse oximetry during full body irradiation, keep the sensor out of the radiation field. If sensor is exposed to the radiation, the reading might be inaccurate or not provided for the duration of the active radiation period.
- Do not use the sensor during MRI scanning or in a MRI environment.
- High ambient light sources such as surgical lights (especially those with a xenon light source), bilirubin lamps, fluorescent lights, infrared heating lamps, and direct sunlight can interfere with the performance of the sensor.
- To prevent interference from ambient light, ensure that the sensor is properly applied, and cover the sensor site with opaque material, if required. Failure to take this precaution in high ambient light conditions may result in inaccurate measurements.
- Elevated levels of Carboxyhemoglobin (COHb) may lead to inaccurate SpO₂ measurements.
- Elevated levels of Methemoglobin (MetHb) will lead to inaccurate SpO₂ measurements.
- Elevated Total Bilirubin levels may lead to inaccurate SpO₂ measurements.
- High levels of COHb or MetHb may occur with a seemingly normal SpO₂. When elevated levels of COHb or MetHb are suspected, laboratory analysis (CO-Oximetry) of a blood sample should be performed.
- Abnormal fingers, intravascular dyes such as indocyanine green or methylene blue or externally applied coloring and texture such as nail polish, acrylic nails, glitter, etc. may lead to inaccurate SpO₂ measurements.
- Inaccurate SpO₂ readings may be caused by severe anemia, low arterial perfusion, motion artifact, hypocapnic or hypercapnic conditions, EMI radiation interference.
- When using multiple Radius PPG sensors, repeat pairing before monitoring to ensure proper wireless connection.
- To prevent damage, do not soak or immerse the sensor in any liquid solution.
- Do not modify or alter the sensor in any way. Alteration or modification may affect performance and/or accuracy.
- Do not attempt to reuse on multiple patients, reprocess, recondition or recycle Masimo sensors or patient cables as these processes may damage the electrical components, potentially leading to patient harm.
- High oxygen concentrations may predispose a premature infant to retinopathy. Therefore, the upper alarm limit for the oxygen saturation must be carefully selected in accordance with accepted clinical standards.
- **Caution:** Replace the sensor when a replace sensor message is displayed, or when a low SIQ message is consistently displayed after completing the low SIQ troubleshooting steps identified in the monitoring device operator's manual.
- **Note:** The sensor is provided with X-Cal® technology to minimize the risk of inaccurate readings and unanticipated loss of patient monitoring. After single-patient use, discard sensor.

INSTRUCTIONS

a) Verification of components

1. Gather all components required for monitoring. Refer to **Fig. 1**.

A. Radius PPG adhesive sensor **B.** Radius PPG wireless bluetooth receiver **C.** Radius PPG reusable chip **D.** Radius PPG reusable chip holder for device

b) Sensor site selection

- Radius PPG ADT: Adult Sensor > 30 kg: The preferred site is the middle or ring finger of the non-dominant hand.
- Radius PPG PDT: Pediatric Sensor 10–50 kg: The preferred site is middle or ring finger of the non-dominant hand.
- Radius PPG INF: Infant Sensor 3–10 kg: The preferred site is the great toe. Alternatively, the toe next to the great toe, or the thumb can be used. 10–20 kg: The preferred site is the middle or ring finger of the non-dominant hand.
- Radius PPG NEO: Neonatal/Adult Sensor < 3 kg: The preferred site is the foot. Alternatively, across the palm and back of the hand can be used. > 40 kg: The preferred site is the middle or ring finger of the non-dominant hand.

c) Attaching the strap to the patient

1. Open the pouch and remove the single-patient-use sensor.
2. Peel off the plastic tab to activate the battery, and discard the tab. Refer to **Fig. 2**.
3. Thread the attachment strap through the plastic loop hole. Wrap the attachment strap around the wrist or limb of the patient (based on selected sensor site) and attach using the hook and loop closure. Refer to **Figs. 3a** and **3b**.
4. Remove the backing from the sensor, if present.

d) Attaching the sensor to the patient.

Adt sensor for ADULTS (> 30 kg) and Pdt sensor for PEDIATRICS (10–50 kg)

1. Orient the sensor so that the detector can be placed first. Place the tip of the finger on the dashed line with the fleshy part of the finger covering the finger outline and detector window. Refer to **Fig. 4a**.
2. Press the adhesive wings, one at a time, onto the finger. Fold the sensor over the finger with the emitter window (*) positioned over the fingernail. Secure the wings down, one at a time, around the finger. Complete coverage of the detector window is needed to ensure accurate data. Refer to **Fig. 4b**.
3. When properly applied, the emitter and detector should be vertically aligned (the black lines should align). Reposition if necessary. Refer to **Fig. 4c**.

Inf sensor for INFANTS (3–10 kg)

1. Direct the sensor cable so that it runs along the top of the foot. Position the detector on the fleshy pad of the great toe. Alternatively, the toe next to the great toe, or the thumb can be used (not shown). Refer to **Fig. 5a**.
2. Wrap the adhesive wrap around the toe so the emitter is positioned on the nail bed of the great toe. Complete coverage of the detector window is needed to ensure accurate data. Refer to **Fig. 5b**.
3. Ensure that the emitter window (*) aligns on the top of the toe directly opposite the detector. Verify correct positioning and reposition if necessary. Refer to **Fig. 5c**.

Neo sensor for NEONATES (< 3 kg)

1. For fragile skin, the stickiness of the medical grade adhesive can be diminished or eliminated by daubing the adhesive areas with a cotton ball or gauze. Refer to **Fig. 6a**.
2. Direct the sensor cable toward the ankle (or wrist). Apply the sensor around the lateral aspect of the foot (or hand), aligned with the fourth toe (or finger). Complete coverage of the detector window is needed to ensure accurate data. Refer to **Fig. 6b**.
3. Wrap the adhesive/foam wrap around the lateral aspect of the foot (or hand) and ensure that the emitter window (*) aligns directly opposite of the detector. Be careful to maintain proper alignment of the detector and emitter windows while attaching adhesive/foam wrap to secure the sensor. Verify correct positioning and reposition if necessary. Refer to **Fig. 6c**.

Neo sensor for ADULTS (> 40 kg) Inf sensor for INFANTS (10–20 kg)

1. Direct the sensor cable so that it runs along the top of the hand. Position the detector on the fleshy part of the finger. Refer to **Fig. 7a**.
2. Alternatively, the sensor may also be applied to the toe (not shown).
3. Wrap the adhesive wrap around the finger so the emitter window (*) aligns on the top of the finger directly opposite the detector. Complete coverage of the detector window is needed to ensure accurate data. Refer to **Fig. 7b**.
4. Check the sensor to verify correct positioning and reposition if necessary. Refer to **Fig. 7c**.

e) Adjusting the sensor cable

1. With the sensor attached on the selected application site, adjust the thin flexible sensor cable. Refer to **Figs. 8a** and **8b**.

f) Reattachment

1. The sensor may be reapplied to the same patient if the emitter and detector windows are clear and the adhesive still adheres to the skin.

g) Applying replacement tape (for Neo sensor only)

1. Remove the existing tape and discard. Refer to **Fig. 9**.
2. Remove the replacement tape from the release liner. Refer to **Fig. 10**.
3. Position the replacement tape over the sensor, aligning the emitter component with the sensor cable. Refer to **Fig. 11**.

h) Disconnecting

1. Push down on the tab to release the reusable chip from the sensor. Refer to **Fig. 12**.
2. After cleaning, store the reusable chip in the chip holder attached to the patient monitor.
3. Discard the adhesive sensor and strap.

SPECIFICATIONS

When used with Masimo SET pulse oximetry monitors, or with licensed Masimo SET pulse oximetry modules the Radius PPG sensors have the following specifications:

Radius PPG sensor used with Masimo device	Radius PPG Adt	Radius PPG Pd1	Radius PPG Inf		Radius PPG Neo	
 Body Weight	> 30 kg	10–50 kg	3–10 kg	10–20 kg	< 3 kg	> 40 kg
Application Site	Finger or Toe	Finger or Toe	Thumb or Great Toe	Finger or Toe	Hand or Foot	Finger or Toe
SpO ₂ Accuracy, No Motion (70–100%) ¹	2%	2%	2%	2%	3% ⁵	2%
SpO ₂ Accuracy, Motion ²	3%	3%	3%	3%	3%	3%
SpO ₂ Accuracy, Low Perfusion ³	2%	2%	2%	2%	3% ⁵	2%
Pulse Rate ¹ Accuracy, No Motion (25–240 bpm)	3 bpm	3 bpm	3 bpm	3 bpm	3 bpm	3 bpm
Pulse Rate Accuracy, Motion ²	5 bpm	5 bpm	5 bpm	5 bpm	5 bpm	5 bpm
Pulse Rate Accuracy, Low Perfusion ⁴	3 bpm	3 bpm	3 bpm	3 bpm	3 bpm	3 bpm

SpO ₂ Upper and Lower Limits of Agreement (LoA)*	
	No Motion
Upper 95% LoA	1.98%
Lower 95% LoA	-2.91%

NOTE: A_{MS} accuracy is a statistical calculation of the difference between device measurements and reference measurements. Approximately two-thirds of the device measurements fell within $\pm A_{MS}$ of the reference measurements in a controlled study.

¹ The Masimo SET Technology has been validated for no motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark pigmented skin in induced hypoxia studies in the range of 70–100% SpO₂ against a laboratory co-oximeter.

² The Masimo SET Technology has been validated for motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark pigmented skin in induced hypoxia studies while performing rubbing and tapping motions, at 2 to 4 Hz at an amplitude of 1 to 2 cm and a non-repetitive motion between 1 to 5 Hz at an amplitude of 2 to 3 cm in induced hypoxia studies in the range of 70–100% SpO₂ against a laboratory co-oximeter.

³ The Masimo SET Technology has been validated for low perfusion accuracy in bench top testing against a Biotek Index 2 simulator and Masimo's simulator with signal strengths of greater than 0.02% and transmission of greater than 5% for saturations ranging from 70% to 100%.

⁴ The Masimo SET Technology has been validated for pulse rate accuracy for the range of 25–240 bpm in bench top testing against a Biotek Index 2 simulator and Masimo's simulator with signal strengths of greater than 0.02% and transmission of greater than 5% for saturations ranging from 70% to 100%.

⁵ The saturation accuracy of the Neonate and Preterm sensors were validated on adult volunteers and 1% was added to account for the properties of fetal hemoglobin.

* See Bland and Altman. Agreement between methods of measurement with multiple observations per individual Journal of Biopharmaceutical Statistics (2007) vol. 17 pp. 571–582.

ENVIRONMENTAL

Storage/Transport Temperature	0–50°C @ ambient humidity
Operating Temperature	0–40°C @ ambient humidity
Storage/Transport Humidity	5–95% non-condensing
Operating Humidity	5–95% non-condensing

BATTERY LIFE

Battery Life	96 hours in typical continuous usage
--------------	--------------------------------------

COMPATIBILITY

 This sensor is intended for use only with devices containing Masimo SET oximetry or pulse oximetry monitors licensed to use Masimo sensors. Radius PPG sensors are compatible with the following Masimo technology boards and software versions and higher: MX version 7.14.8.x. Each sensor is designed to operate correctly only on the pulse oximetry systems from the original device manufacturer. Use of this sensor with other devices may result in no or improper performance.

For Compatibility Information Reference: www.Masimo.com

WARRANTY

Masimo warrants to the initial buyer only that these products, when used in accordance with the directions provided with the Products by Masimo, will be free of defects in materials and workmanship for a period of six (6) months. Single use products are warranted for single patient use only.

THE FOREGOING IS THE SOLE AND EXCLUSIVE WARRANTY APPLICABLE TO THE PRODUCTS SOLD BY MASIMO TO BUYER. MASIMO EXPRESSLY DISCLAIMS ALL OTHER ORAL, EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. MASIMO'S SOLE OBLIGATION AND BUYER'S EXCLUSIVE REMEDY FOR BREACH OF ANY WARRANTY SHALL BE, AT MASIMO'S OPTION, TO REPAIR OR REPLACE THE PRODUCT.

WARRANTY EXCLUSIONS

This warranty does not extend to any product that has been used in violation of the operating instructions supplied with the product, or has been subject to misuse, neglect, accident or externally created damage. This warranty does not extend to any product that has been connected to any unintended device or system, has been modified, or has been disassembled or reassembled. This warranty does not extend to sensors or patient cables that have been reprocessed, conditioned or recycled.

IN NO EVENT SHALL MASIMO BE LIABLE TO BUYER OR ANY OTHER PERSON FOR ANY INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING WITHOUT LIMITATION LOST PROFITS), EVEN IF ADVISED OF THE POSSIBILITY THEREOF. IN NO EVENT SHALL MASIMO'S LIABILITY ARISING FROM ANY PRODUCTS SOLD TO BUYER (UNDER A CONTRACT, WARRANTY, TORT OR OTHER CLAIM) EXCEED THE AMOUNT PAID BY BUYER FOR THE LOT OF PRODUCT(S) INVOLVED IN SUCH CLAIM. IN NO EVENT SHALL MASIMO BE LIABLE FOR ANY DAMAGES ASSOCIATED WITH A PRODUCT THAT HAS BEEN REPROCESSED, RECONDITIONED OR RECYCLED. THE LIMITATIONS IN THIS SECTION SHALL NOT BE DEEMED TO PRECLUDE ANY LIABILITY THAT, UNDER APPLICABLE PRODUCTS LIABILITY LAW, CANNOT LEGALLY BE PRECLUDED BY CONTRACT.








NO IMPLIED LICENSE

THIS SINGLE-PATIENT SENSOR IS LICENSED TO YOU UNDER THE PATENTS OWNED BY MASIMO FOR SINGLE-PATIENT USE ONLY. BY ACCEPTANCE OR USE OF THIS PRODUCT, YOU ACKNOWLEDGE AND AGREE THAT NO LICENSE IS GRANTED FOR USE OF THIS PRODUCT WITH MORE THAN A SINGLE PATIENT.

AFTER SINGLE-PATIENT USE, DISCARD SENSOR. PURCHASE OR POSSESSION OF THIS SENSOR CONFERS NO EXPRESS OR IMPLIED LICENSE TO USE THE SENSOR WITH ANY SINGLE WHICH IS NOT SEPARATELY AUTHORIZED TO USE MASIMO SENSORS.

CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events. The following symbols may appear on the product or product labeling:

SYMBOL	DEFINITION	SYMBOL	DEFINITION	SYMBOL	DEFINITION
	Follow instructions for use		Separate collection for electrical and electronic equipment (WEEE)	Rx ONLY	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician
	Consult instructions for use		Lot code		Mark of conformity to European Medical Device Directive 93/42/EEC
	Manufacturer		Catalogue number (model number)		Authorized representative in the European community
	Date of manufacture YYYY-MM-DD		Masimo reference number		Body weight
	Use-by YYYY-MM-DD		Light Emitting Diode (LED) LED emits light when current flows through		Storage temperature range
	Fragile, handle with care		Greater than		Keep dry
	Non-Sterile		Less than		Do not use if package is damaged
	Not made with natural rubber latex		Storage humidity limitation		Atmospheric pressure limitation
	Caution		Bluetooth	IP22	Protected against solid foreign objects of 12.5 mm diameter and greater and protection against vertically falling water drops when enclosure is tilted at 15 degrees
	Do not re-use/Single patient use only		Federal Communications Commission (FCC) Licensing	FCC ID:	Identifies unit has been registered as a radio device
	Instructions/Directions for Use/Manuals are available in electronic format @ http://www.Masimo.com/TechDocs Note: eIFU is not available in all countries.				

Patents: <http://www.masimo.com/patents.htm>

Radius PPG is a trademark of Masimo Corporation.

Masimo, SET, X-Cal, and  are federally registered trademarks of Masimo Corporation.

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PERFORMANCE SPECIFICATIONS


The table below shows A_{RMS} (Accuracy Root Mean Square) values measured using the Radius PPG sensors under no motion, with Masimo Technology in a clinical study.

Radius PPG Sensor	
SpO ₂	A_{RMS}
90-100%	1.14 (1,14) %
80-90%	1.29 (1,29) %
70-80%	1.41 (1,41) %
70-100%	1.33 (1,33) %

SpO ₂ Upper and Lower Limits of Agreement (LoA)*	
	Actual Value
Upper 95% LoA	1.98%
Lower 95% LoA	-2.91%

*See Bland and Altman. Agreement between methods of measurement with multiple observations per individual *Journal of Biopharmaceutical Statistics* (2007) vol. 17 pp. 571-582.

70-100%
 (SpO₂-SaO₂) vs. (SpO₂+SaO₂)/2 Bland Altman fit and upper 95% and lower 95% limits of agreement.

© 2019 Masimo Corporation

 **Manufacturer:**
Masimo Corporation
52 Discovery
Irvine, CA 92618
USA

www.masimo.com

EU Authorized Representative for Masimo Corporation:

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MDSS GmbH